

SHUKLATM

Medical

Innovators of
Xtract-All[®]

Universal Broken, Stripped and Seized Screw Extraction System

Removes:

- M2 to M8 Size Broken & Stripped Screws
 - Cold Weld (fused) Stripped Screws
- Capable of applying over 250 in-lbs of Torque*

Attachments include:



- 48 Extractors
- 3 Extensions
- Breaker Bar
- Jacobs Chuck Attachment
- Ratcheting T Handle



*Broken Screw
Extraction System*



*Stripped and
Seized Screw
Extraction System*




**Have the instrument you need
when you need it.**

Contact us at: **888-4shukla (888-474-8552)**
or **customer.service@ShuklaMedical.com**
See our website at: **www.ShuklaMedical.com**



Universal Broken, Stripped and Seized Screw Extraction System

SHUKLATM
Medical



Kit: MKS 1022, Case: CA145, Lid: CA146

(MKS 1022) BROKEN, STRIPPED & SEIZED BONE SCREW EXTRACTION SYSTEM

No.	Part Number	Description	Units	Single Use
1	TP341-M01	Extractor tip for 2.0-2.5 mm Broken Screw	4	✓
2	TP342-M01	Extractor tip for 2.5-4.0 mm Broken Screw	4	✓
3	TP343-M01	Extractor tip for 3.5-5.0 mm Broken Screw	4	✓
4	TP344-M01	Extractor tip for 4.5-6.0 mm Broken Screw	4	✓
5	TP345-M01	Extractor tip for 5.5-7.0 mm Broken Screw	4	✓
6	TP346-M01	Extractor tip for 7.0-8.0 mm Broken Screw	4	✓
7	TP347-M01	Extractor tip for 2.0-2.5 mm Stripped Screw Head Recess	4	✓
8	TP348-M01	Extractor tip for 2.5-3.5 mm Stripped Screw Head Recess	4	✓
9	TP349-M01	Extractor tip for 3.5-4.5 mm Stripped Screw Head Recess	4	✓
10	TP350-M01	Extractor tip for 4.5-5.5 mm Stripped Screw Head Recess	4	✓
11	TP351-M01	Extractor tip for 5.5-6.5 mm Stripped Screw Head Recess	4	✓
12	TP352-M01	Extractor tip for 6.0-7.0 mm Stripped Screw Head Recess	4	✓
13	MMI1152-M01	4" Extension Driver	1	
14	MMI1153-M01	6" Extension Driver	1	
15	MMI1154-M01	8" Extension Driver	1	
16	CA145	Case, Stripped and Seized Screw Removal System	1	
17	CA146	Lid, Stripped and Seized Screw Removal System	1	
18	CA139	Caddy, Stripped and Seized Screw Removal System	1	
19	CA140	Caddy lid, Stripped and Seized Screw Removal System	1	
20	CA143	Caddy, Broken Screw Extraction System	1	
21	CA144	Caddy lid, Broken Screw Extractor System	1	
22	HD233-M01	Asymmetric Cannulated Ratcheting T-Handle	1	
23	AD120-M01	Power Drill Adapter	1	
24	HD239-M01	Breaker Bar	1	
25	HD239-M01	Stripped Screw Extractor Breaker Bar	1	

