

Xtract-All[®] mini

Universal Small Screw Extraction System

Designed to remove small bone screws commonly used in hand, foot, cranio-maxillo-facial and spine procedures

...Because Every Minute Counts[™]




- ✦ **Wide Variety of Drivers for Efficient Removals**
- ⬡ **Removes Broken, Stripped and Seized Bone Screws**
- ✦ **Ratcheting Handle Designed for Enhanced Driver Control**

Xtract-All Mini® (S9MINI) Universal Ratcheting Screwdriver Set

No.	Part Number	Description	Qty
1	SCS201	Case	1
2	SHN201	Ratcheting Handle	1
3	SDR201	Hex 1.5mm	1
4	SDR202	Hex 2.0mm	1
5	SDR203	Hex 2.5mm	1
6	SDR204	Small Cruciform	1
7	SDR205	Medium Cruciform	1
8	SDR206	Large Cruciform	1
9	SDR207	Slotted 1.5mm - 2.0mm	1
10	SDR208	Slotted 3.0mm	1
11	SDR209	Center 1.0mm	1
12	SDR210	Center 1.5mm	1
13	SDR211	Center 2.0mm	1
14	SDR212	Center 2.7mm	1

No.	Part Number	Description	Qty
15	SDR213	Center 3.5mm	1
16	SDR214	Clutch 2.7mm	1
17	SDR215	Phillips 2.0mm - 2.7mm	1
18	SDR216	T4 Torx	1
19	SDR217	T5 Torx	1
20	SDR218	T6 Torx	1
21	SDR219	T7 Torx	1
22	SDR220	T8 Torx	1
23	SDR221	T10 Torx	1
24	SDR222	Osteotome	1
*25	SXT201	Broken Screw Extractor	2
*26	SXT202	Small Stripped Screw Extractor	2
*27	SXT203	Large Stripped Screw Extractor	2



SXT201, SXT202, and SXT203 are single use only and must be discarded after use.

For additional information and pricing, please contact customer service.

Phone & Fax: 732-474-1769 | Email: Customer.Service@ShuklaMedical.com

Website: www.ShuklaMedical.com | Address: 151 Old New Brunswick Road, Piscataway NJ 08854