

WINQUIST 4

UPGRADE Xtract-Air[®] NAIL

Universal IM Nail Extraction System

The first and most comprehensive IM nail extraction system on the market.
Now in its 4th generation.



58% INCREASE IN EXTRACTION INSTRUMENTATION



EXPANDED RANGE OF EXTRACTOR SIZES make threaded tips, solid IM nail extractors, & conical extractors more effective.

[See reverse for more >](#)



SHUKLA MEDICAL[®]
Universal Orthopedic Extraction Technologies

...Because Every Minute Counts[®]


BONE-CLEARING POWER
BONE-CLEARING FLUTES on threaded extractor tips and **TREPINE TEETH** on solid IM nail extractors help ease attachment



NO SECOND-GUESSING
THREAD SIZING BOARD matches a threaded IM nail cap to a Winquist threaded extractor



ERGONOMIC ADJUSTMENTS
MODULAR C-FRAME adjusts to user preference with an optional, adjustable over-handed hold



RATCHETING WRENCH makes tightening the hex nut on the c-frame easier than ever



SILICONE COATED MALLET absorbs shock and can be swung with two hands for maximum controlled impaction force



WINQUIST 4


UPGRADE Xtract-All[®] NAIL

Universal IM Nail Extraction System


58% INCREASE IN EXTRACTION INSTRUMENTATION

28 THREADED TIPS • 5 SOLID EXTRACTORS • 3 CONICAL EXTRACTORS • 2 HOOK EXTRACTORS
bringing the total extractor count from 24 (Winquist 3) to 38 (Winquist 4)

Winquist 3 (MKS1017)



Winquist 4 (S9WINQ)



Previous Winquist 3 (MKS1017) versus Upgraded Winquist 4 (S9WINQ)

Winquist 3 MKS1017	Winquist 4 S9WINQ	Part Description
	✓	8-32 Threaded Extractor
	✓	10-24 Threaded Extractor
	✓	10-32 Threaded Extractor
	✓	12-28 Threaded Extractor
	✓	12-24 Threaded Extractor
✓	✓	1/4-20 Threaded Extractor
✓	✓	1/4-28 Threaded Extractor
✓	✓	5/16-18 Threaded Extractor
✓	✓	5/16-24 Threaded Extractor
	✓	5/16-32 Threaded Extractor
✓	✓	3/8-16 Threaded Extractor
✓	✓	3/8-24 Threaded Extractor
	✓	7/16-14 Threaded Extractor
✓	✓	7/16-20 Threaded Extractor
✓	✓	M5 x 0.8 Threaded Extractor
✓	✓	M5.5 x 0.9 Threaded Extractor
✓	✓	M6 x 0.75 Threaded Extractor
✓	✓	M6 x 1.0 Threaded Extractor
	✓	M7 x 1.0 Threaded Extractor
✓	✓	M8 x 1.0 Threaded Extractor
✓	✓	M8 x 1.25 Threaded Extractor
	✓	M10 x 1.0 Threaded Extractor
	✓	M10 x 1.25 Threaded Extractor
✓	✓	M10 x 1.5 Threaded Extractor
	✓	M12 x 1.25 Threaded Extractor
✓	✓	8-32 Female Threaded Extractor
✓	✓	1/4-20 Female Threaded Extractor
✓	✓	3/8-32 Female Threaded Extractor
one-size-fits-all conical extractor	✓	Small Conical Extractor (4.0mm-8.0mm)
one-size-fits-all conical extractor	✓	Medium Conical Extractor (8.0mm-12mm)
one-size-fits-all conical extractor	✓	Large Conical Extractor (12mm-16mm)
✓	✓	Solid IM Nail Extractor (6.0mm-8.0mm)
✓	✓	Solid IM Nail Extractor (8.0mm-10mm)
✓	✓	Solid IM Nail Extractor (10mm-12.5mm)
✓	✓	Solid IM Nail Extractor (12mm-15mm)
	✓	Solid IM Nail Extractor (15mm-18mm)
✓	✓	Small Hook
✓	✓	Large Hook



...Because Every Minute Counts[®]